




Con@erto

*The combined Mesotherapy, PRP
& Carboxytherapy*



Mésotherapy, PRP & Carboxytherapy : All in one !

Synthesis of all existing injection techniques, the Concerto is the result of more than 20 years of experience in the field of medical injectors.

Perfect combination of Mesotherapy and Carboxytherapy, the Concerto manages programs face, body and scalp.

Intuitive, light and accurate it injects the right dose, in the right place and at the right depth. All types of drug and aesthetic solutions can be used and also the carbon dioxide gas.



Efficient and easy to use

- **PRECISE** : Electronic control of the depth, angle, speed and volume of each injection.
- **INTUITIVE** : An exclusive system adjusts automatically the gas flow or quantity of product to the kind of the needle used.
- **AUTOADAPTATIVE** : Automatic adjustment of the volume delivered according to the injection speed selected.
- **PAINLESS** : Lightness electronic needle support for quick and comfortable injections in contact with the skin.
- **COMPLETE** : 2 injection modes
 - **Manual** for the face, eyes, wrinkles...
 - **Automatic** for the body, abdomen, legs...
- **UNIVERSAL** : All Mesotherapy and PRP products ; all kind of CO2 sources: Bottle or cartridge.

Mesotherapy and PRP



Mesotherapy treats in situ by injection through a fine needle. Suitable for precise treatments, it offers the possibility of a localized application with the exact dose of active products.

Without global diffusion inside the organism the Mesotherapy injection allows the use of concentrated active solution and the treatment of pathologies like rheumatology, arthrosis, tendinitis, sport medicine and aesthetic medicine.

Mesotherapy for aesthetic applications

For aesthetic treatments, Mesotherapy is intended for injections of hyaluronic acid, vitamins, PRP or fillers.

The main indications are wrinkles, hair treatments, increase of volume (lips, cheeks), stretch of the skin, anti-aging.



3 COMPLEMENTARY PROGRAMS FOR MESOTHERAPY

- Nappage :** For precise and regular lines of injections. Automatic adaptation of the dose delivered according to the injection speed selected and the needle used.
- Expert :** Special program dedicated for point per point and PRP injections. Adjustable depth (from 1 to 10 mm) and dose (from 0.01 ml to infinity). Optimization of the quantity injected and the protocols effectiveness.
- Fillers :** Perfectly regular retro-injections suitable for fillers and volume.


Carboxytherapy

The Carboxytherapy® is the subcutaneous injection of carbon dioxide (CO₂ gas). Naturally generated by the organism as a metabolism consequence, the CO₂ increases immediately the flow of blood and oxygen inside the muscular and adipose tissues.

Locally injected the carbonic gas is automatically removed by the blood circulation and eliminated by the lungs. In front of CO₂, the hemoglobin balances the O₂ and CO₂ levels by release of oxygen (Bohr effect).

Performed several times, the Carboxytherapy increases the microcirculation and tissues oxygenation, it breaks the fibrosis and restructure the adipose tissues for therapeutic, rheumatological and aesthetic treatments with natural and lasting results.

Carboxytherapy for aesthetic applications



For aesthetic treatments the Carboxytherapy improves skin elasticity and tissues vitality. Bring of oxygen done by Carboxytherapy improves almost all anti-aging possibilities in a natural and harmonious way.

Carboxytherapy for aesthetic is mainly used for the treatment of wrinkles, stretch marks, scars, cutaneous relaxing (décolleté, back of the arms...), cellulite (double chin, oval of the face, thighs, abdomen, knees...), hair loss and dark circles.



Suitable for all Carboxytherapy treatments



Automatic : For body treatments such as stretch marks or skin relaxing.

Manual : For injections on sensitive areas such as eyes area and mouth.

Flux regulation : Patented adjustment of gas flow according to the needle used and the area treated.

General features

Total autonomy	Works indifferently on battery with 20 hours of autonomy (Li-ion) or connected to the main supply (90V - 240V).
Friendly use	High-definition graphics screen, touch pad control, ready for session in a few seconds.
Ergonomic	Exclusive design that gives to the Concerto® lightness and easy to use application. Ambidextrous use.
Painless	Needle projection speed smoothly, no phenomenon painful vibration induced.

Briefly

- Directional **lighting (LED)** of the treated area.
- Simultaneous or post-injection micro massage for better products diffusion.
- **Several kind of stabilizers** available for all treatments areas.
- No exclusive kit, the Concerto® can be used with **all needles, products and syringes**.
- **Automatic detection** of the syringe and viscosity of products injected.
- Designed and manufactured in France.